

Setting Up Your Messages

The iGlucose® Virtual Coach sends tips, reminders and trends via text or email to your selected method of contact.

1. If you would like text messages, please enter your phone number and the confirmation code sent to your phone.

Add Mobile/Cell Phone Number

+1 xxx-xxx-xxxx

Enter Confirmation Code

2. If you would like email, please enter your email address and the confirmation code sent to your email.

Add Email

jane.doe@yahoo.com

Enter Confirmation Code


3. Choose the type of messages you would like to receive.

- | | |
|--|--|
| <input type="checkbox"/> BG Result Message | <input type="checkbox"/> Critical Range Messages |
| <input type="checkbox"/> Daily Recap | <input type="checkbox"/> Tips and Reminders |

Create Your Circle of Care

Invite your friends, family, doctors and diabetes nurses to see how you are doing and to help you when you need it.

They can choose if they want to receive texts, emails or just access to your portal.



1. Select if the member of your Circle of Care is a parent, child, healthcare provider, etc.

Who is the recipient?

Please Select

- Parent of Person with Diabetes
- Child of Person with Diabetes
- Other Caregiver (family member or friend)
- Doctor
- Diabetes Educator
- Insurer

2. Enter their email address so we can invite them to receive messages when the meter is used.

Be sure to tell them to look out for a confirmation email from support@iglucose.com to activate their account.

Enter their email so we can invite them to receive messages:

jane.doe@yahoo.com

We Hope You Enjoy Using the iGlucose® Cellular Diabetes Care Solution

For assistance, please contact us at:

1-844-iGlucose (1-844-445-8267)

or Email us at support@iglucose.com



iGlucose®
Blood Glucose Monitoring System

iGlucose®

Blood Glucose Monitoring System

Quick Start Guide

Website Portal, Texting
and Circle of Care

Introduction

Welcome to the iGlucose®
Cellular Diabetes Care Solution

The iGlucose System was designed so that you can get the support you need in managing your diabetes. By connecting your meter to the portal, you and your designated caregivers can receive updates on your progress via email and/or text.

You and your caregivers can also log on to the web portal to view logbooks and graphs of your blood glucose readings.


Let's Get Started Now! ➤

Setting Up Your Account

Setting up your account and connecting your iGlucose® System allows you to receive tips and reminders from the iGlucose Coach, login to the website portal to see your charts and trends, and for your caregivers to receive messages on how you are doing.

Connecting Your Meter

- 1. Go to connect.igluco.se.com.
- 2. Enter your Device ID (found on the label located on the back of the meter).



Enter iGlucose Device ID:
3088470

3. Enter the name of the person with diabetes.

| | |
|------------|-----------|
| First Name | Last Name |
| Jane | Doe |

4. Enter your time zone.
So iGlucose® Coach knows when to send you reminders

Please Select Time Zone

5. Enter a username and password.

| | | |
|-------------------|-----------------|-------------------|
| Email or Username | Create Password | Re-enter Password |
| jdoe | ocean | ocean |

6. Personal information and method of receiving messages.

| | |
|--------------------|-----------------|
| Email | Phone Number |
| jane.doe@yahoo.com | +1 xxx-xxx-xxxx |

The iGlucose Coach would like to help you manage your diabetes by sending text or email messages.

Which delivery method would you prefer for your messages?

Please Select Method

Text
Email

Complete Your Diabetes-Related Information

1. If your doctor has recommended a blood glucose range, enter it.

Blood Glucose Range

| | | |
|-----|------|-----------------|
| Low | High | After Meal High |
| 80 | 130 | 180 |

Note: The default range is based on the American Association of Diabetes.

2. Enter your recommended test frequency.

Tests per Day

2

3. Set up reminders, so you don't miss checking your BG.

Note: Reminders will be sent at the specified time in the time zone previously identified above.

Testing Times

| | | | |
|---------|---|---------|---|
| 8:00 am | 🕒 | 8:00 pm | 🕒 |
|---------|---|---------|---|

☒ Send me reminders to test

