


iDigiTemp

Cellular-Connected Infrared Thermometer

The First Handheld Cellular Thermometer for RPM

Revolutionizing how temperature data is captured and transmitted in real time. iDigiTemp seamlessly integrates with Smart Meter's Private Data Network, ensuring secure, reliable RPM connectivity and enabling instant transmission of accurate temperature readings to your RPM platform.

- A Remote Patient Monitoring (RPM) solution with real-time testing results, securely transmitted to the care team.
- Clinicians receive real-time data about patient temperature, a vital sign that often signals early changes in a patient's condition.
- Enhanced cellular connectivity with new multicarrier SIM technology to ensure vital sign data is delivered securely and consistently.
- Bright, easy-to read LCD screen.



Connected by



Features

- FDA-Cleared
- Enhanced cellular connectivity to our Private Data Network via multicarrier SIM technology, connecting to the strongest cellular network available in the device area
- Temperature unit Celsius (C) or Fahrenheit (F)
- Measurements from forehead or ear site
- Indicator color and sound alerts with mute/un-mute capability
- Advanced infrared temperature sensor with high sensitivity
- Automatic switching from Forehead mode to Ear mode when taking the probe cap on and off
- Memory Stores up to 20 measurements
- One-button operation
- One-second measurement time
- Powered by two AAA batteries
- Automatic power off in 10 seconds

Why Smart Meter is the No. 1 Supplier of Cellular RPM™ Solutions

The Only Trusted Private Data Network for RPM in the U.S.

- Smart Meter's Private Data Network for RPM is an exclusive ecosystem with robust security layers protecting patient-generated health data from our proprietary monitoring devices to our SmartPartners™ platforms.



Cellular RPM™ is Better Than Bluetooth®

- Patients have problems using Bluetooth – it's that simple.
- Our proprietary cellular-enabled devices deliver the most patient-friendly experience available supporting more frequent testing, additional revenue, an improved volume of health data facilitating care management, and higher customer and patient satisfaction.



70% of people have trouble sending health data with BlueTooth® devices.

Everything Working Seamlessly in the Background

- We simplify the deployment and management of remote patient monitoring so you can focus on transforming healthcare.
- SmartPartners have exclusive access to curated monitoring devices, real-time RPM data stacks, shipping and logistics management, automated refills, and success training services.



100% Quality Guarantee

- We stand behind our solutions: products, data, platform, and services.
- If anything fails to meet your standards, we will resolve the problem to your satisfaction. Guaranteed!



Grow With Us

Whether you're implementing an RPM system for the first time or need to reach new customers with expanded solutions, we're here to help.

www.smartmeterrpm.com | 1-813-773-4080



SMARTMETER™
Engaging Health Data

6005 Benjamin Rd
Suite A
Tampa, FL 33634

DT-S-IDTPGR-U-1025