

# iGlucose® Plus


## For Management of Gestational Diabetes

Gestational diabetes can develop during pregnancy in women who don't already have diabetes. Every year, 5% to 9% of U.S. pregnancies are affected by gestational diabetes.\*

## The Gold Standard in Cellular Blood Glucose Monitoring

The new and improved iGlucose® Plus sets the standard for accuracy, reliability, and convenience, in Cellular Blood Glucose Monitoring for diabetes management.

- A Remote Patient Monitoring (RPM) solution with real-time testing results, securely transmitted to the care team.
- Gold standard in blood glucose testing includes real gold in the test strips for more accurate results.
- Enhanced cellular connectivity with new multicarrier SIM technology for reliable transmittal of vital signs to care team.
- Bright, easy-to read colored LCD screen, with the option to view in English or Spanish.
- Customizable on-screen alerts and reminders to help patients test more.
- Approved for use during pregnancy\*\*



Connected by



## Features

- FDA-Cleared
- Blood Glucose Range of 20-600 mg/dL
- Easy to Read, 36.6 mm X 49 mm Large Colored LCD Display Screen (English or Spanish)
- Enhanced cellular connectivity to our Private Data Network via multicarrier SIM technology, connecting to the strongest cellular network available in the device area
- Customizable On-Screen Alerts and Reminders (English or Spanish)
- HCT Supports a Range Between 20% and 60%
- Automatic Strip Coding
- Memory Stores up to 500 Blood Glucose Tests with Date and Time
- Meal Markers can be Set On/Off on the Display Screen for Pre/Post-Meal
- 5-Second Test Time
- Rechargeable 3.7V Li-Ion Battery
- Includes Airplane Mode Setting for Easy Use While Traveling
- Result Averages Shown for the Last 7, 14, 30, 60 and 90 Days of Testing

## Pricing

**\$299.95 for 3 months,  
After 3 months, \$15/month**

*Includes iGlucose Plus kit with iGlucose Plus meter, lancing device, lancets, unlimited test strips, shipping directly to patient and access to ordering and data portal*



\* 2024 Centers for Disease Control and Prevention statistics  
\*\* CGM is not currently FDA approved for use during pregnancy in the US



# Why Smart Meter is the No. 1 Supplier of Cellular RPM™ Solutions

## The Only Trusted Private Data Network for RPM in the U.S.

- Smart Meter's Private Data Network for RPM is an exclusive ecosystem with robust security layers protecting patient-generated health data from our proprietary monitoring devices to our SmartPartners™ platforms.



## Cellular RPM™ is Better Than Bluetooth®


- Patients have problems using Bluetooth – it's that simple.
- Our proprietary cellular-enabled devices deliver the most patient-friendly experience available supporting more frequent testing, additional revenue, an improved volume of health data facilitating care management, and higher customer and patient satisfaction.



70% of people have trouble sending health data with BlueTooth® devices.

## Everything Working Seamlessly in the Background

- We simplify the deployment and management of remote patient monitoring so you can focus on transforming healthcare.
- SmartPartners have exclusive access to curated monitoring devices, real-time RPM data stacks, shipping and logistics management, automated refills, and success training services.



## 100% Quality Guarantee

- We stand behind our solutions: products, data, platform, and services.
- If anything fails to meet your standards, we will resolve the problem to your satisfaction. Guaranteed!



## Grow With Us

Whether you're implementing an RPM system for the first time or need to reach new customers with expanded solutions, we're here to help.



**SMARTMETER**  
Engaging Health Data

6005 Benjamin Road  
Suite A  
Tampa, FL 33634

S-S-GESTDIA-U-0325

[www.smartmeterrpm.com](http://www.smartmeterrpm.com) | 1-813-773-4080