


iGlucose[®] Touch

Cellular-Connected Blood Glucose Monitoring System

A Touch-Enabled Glucose Meter Built for RPM

iGlucose Touch brings a new level of simple, intuitive testing with a touchscreen design and built-in cellular connectivity. Every test provides fast and accurate results that are instantly transmitted through Smart Meter's Private Data Network securely to your RPM platform.

- A Remote Patient Monitoring (RPM) solution with real-time testing results, securely transmitted to the care team.
- A bright, responsive touchscreen for easy navigation, delivering a user-friendly experience without the need for multiple buttons.
- Enhanced cellular connectivity with new multicarrier SIM technology to ensure vital sign data is delivered securely and consistently.
- Markers with each test result allow for easy recording of relevant statuses and events.



Connected by



Features

- FDA-Cleared
- Enhanced cellular connectivity to our Private Data Network via multicarrier SIM technology, connecting to the strongest cellular network available in the device area
- Blood Glucose Range of 20-600 mg/dL
- 2.8 in Color LCD touchscreen display
- Automatic strip coding
- Memory stores up to 1000 blood glucose readings with dates and times
- Meal markers can be set with a target blood glucose range and displayed next to measurement value
- 5-Second test time
- Rechargeable lithium battery (3.7V)
- Automatic turn off time can be set by user or default to turn off after 30 seconds

www.smartmeterrpm.com | 1-813-773-4080



SMARTMETER™
Engaging Health Data

Why Smart Meter is the No. 1 Supplier of Cellular RPM™ Solutions

The Only Trusted Private Data Network for RPM in the U.S.

- Smart Meter's Private Data Network for RPM is an exclusive ecosystem with robust security layers protecting patient-generated health data from our proprietary monitoring devices to our SmartPartners™ platforms.



Cellular RPM™ is Better Than Bluetooth®

- Patients have problems using Bluetooth – it's that simple.
- Our proprietary cellular-enabled devices deliver the most patient-friendly experience available supporting more frequent testing, additional revenue, an improved volume of health data facilitating care management, and higher customer and patient satisfaction.



70% of people have trouble sending health data with BlueTooth® devices.

Everything Working Seamlessly in the Background

- We simplify the deployment and management of remote patient monitoring so you can focus on transforming healthcare.
- SmartPartners have exclusive access to curated monitoring devices, real-time RPM data stacks, shipping and logistics management, automated refills, and success training services.



100% Quality Guarantee

- We stand behind our solutions: products, data, platform, and services.
- If anything fails to meet your standards, we will resolve the problem to your satisfaction. Guaranteed!



Grow With Us

Whether you're implementing an RPM system for the first time or need to reach new customers with expanded solutions, we're here to help.

www.smartmeterrpm.com | 1-813-773-4080



SMARTMETER™
Engaging Health Data

6005 Benjamin Rd
Suite A
Tampa, FL 33634

IG-S-IGTPGR-U-1125