

ATHENAHEALTH API INTEGRATION ONBOARDING GUIDE


Pre-Requisites

Business or Care Management are required.
 Signed Athena Addendum (Pricing Agreement).

- ▶ Fees are billed per transaction with Athena API
- ▶ Addendum should include the practice ID (see below on how to obtain)

Obtaining Practice ID

Log in to your Athena account.
 Click on your name, top right corner.
 Copy the Practice ID from the menu (see below).



Getting Started

Customer must log into Athena and be assigned one the following athenaOne roles:

- ▶ Communicator admin
- ▶ Practice Superuser
- ▶ Superuser - Clinicals
- ▶ Superuser - Collector
- ▶ Superuser - Communicator
- ▶ Management - Clinicals
- ▶ Management - Collector
- ▶ Management - Communicator

Navigate to the Manage Certified API Access admin page.

- ▶ On the Main Menu, click the “Settings” icon
- ▶ Under ADMIN, click Practice Manager
- ▶ In the Task Bar, under PRACTICE LINKS — Certified API Apps, click Manage Certified API Access

The Manage Certified API App Access page appears.

Search the table to view all apps that have access to your tablespace. Use the search bar to find “Smart Meter RPM”.

Use the toggle on the right to enable or disable the app (blue = enabled, grey = disabled).

If you use the toggle to enable or disable an app, a pop up will appear to confirm you would like to complete this action, select Continue.